




# Closed Suction System

陕械注准20232080067



# Product features


- The sputum suction tube is completely wrapped in a film cover to avoid pollution to the surrounding environment and medical personnel.
- Manual negative pressure suction device. It is equipped with flush and administration connectors.
- The effective time of 24 hours or 72 hours for use can be selected for different products.
- Individual packaging, ethylene oxide sterilization, single use.



# Clinical guideline

- Used for endotracheal tube for patients undergoing general anesthesia or critical patients for sputum suction use in rescue.
- Used for sputum suction for tracheotomy patients under continuous oxygen supply.


## Closed suction system



DRMXPC01			
Code	Size	O.D	Qty
081704	12F	4mm	500
081705	14F	4.67mm	500
081706	16F	5.33mm	400

\*For more models and specifications, please refer to the product search table.


## Closed suction system



DRMXCCT01			
Code	Size	O.D	Qty
081713	10F	3.33mm	500
081714	12F	4mm	500
081715	14F	4.67mm	500

\*For more models and specifications, please refer to the product search table.


## Closed suction system



DRMXCPC02			
Code	Size	O.D	Qty
081709	12F	4mm	500
081710	14F	4.67mm	500
081711	16F	5.33mm	400

\*For more models and specifications, please refer to the product search table.

## Closed suction system



DRMXCCT02			
Code	Size	O.D	Qty
081718	10F	3.33mm	500
081719	12F	4mm	500
081720	14F	4.67mm	500

\*For more models and specifications, please refer to the product search table.

# Product features


- The sputum suction tube is completely wrapped in a film cover to avoid pollution to the surrounding environment and medical personnel.
- Manual negative pressure suction device. It is equipped with flush and administration connectors.
- The effective time of 24 hours or 72 hours for use can be selected for different products.
- Individual packaging, ethylene oxide sterilization, single use.



# Clinical guideline

- Used for endotracheal tube for patients undergoing general anesthesia or critical patients for sputum suction use in rescue.
- Used for sputum suction for tracheotomy patients under continuous oxygen supply.


## Closed suction system



DRMXCCZ01			
Code	Size	O.D	Qty
081723	10F	3.33mm	500
081724	12F	4mm	500
081725	14F	4.67mm	500


\*For more models and specifications, please refer to the product search table.

## Closed suction system



DRMXPR			
Code	Size	O.D	Qty
081700	6F	2mm	500
081701	8F	2.67mm	500

## Closed suction system



DRMXCCZ03			
Code	Size	O.D	Qty
081728	10F	3.33mm	500
081729	12F	4mm	500
081730	14F	4.67mm	500

\*For more models and specifications, please refer to the product search table.